



Contents lists available at ScienceDirect

Journal of the Mechanical Behavior of Biomedical Materials

journal homepage: www.elsevier.com/locate/jmbbm

Biomechanical testing of zirconium dioxide osteosynthesis system for Le Fort I advancement osteotomy fixation

Lukas Hingsammer^{a,c,*}, Markus Grillenberger^b, Martin Schagerl^b, Malek Michael^a, Hunger Stefan^a^a Kepler University Clinic Linz, Department of Maxillofacial Surgery, Krankenhausstrasse 9, 4020 Linz, Austria^b Johannes Kepler University Linz, Institute of Constructional Lightweight Design, Altenbergstrasse 69, 4040 Linz, Austria^c University Clinic of Zurich, Department of Maxillofacial Surgery, Krankenhausstrasse 24, 8032 Zurich, Switzerland

ARTICLE INFO

Keywords:

Maxillofacial surgery
Patient specific implants
Zirconium dioxide
Osteosynthesis system

ABSTRACT

The following work is the first evaluating the applicability of 3D printed zirconium dioxide ceramic miniplates and screws to stabilize maxillary segments following a Le-Fort I advancement surgery. Conventionally used titanium and individual fabricated zirconium dioxide miniplates were biomechanically tested and compared under an occlusal load of 120 N and 500 N using 3D finite element analysis. The overall model consisted of 295,477 elements. Under an occlusal load of 500 N a safety factor before plastic deformation respectively crack of 2.13 for zirconium dioxide and 4.51 for titanium miniplates has been calculated.

From a biomechanical point of view 3D printed ZrO₂ mini-plates and screws are suggested to constitute an appropriate patient specific and metal-free solution for maxillary stabilization after Le Fort I osteotomy.

1. Introduction

Le Fort I osteotomy is a well-established surgical technique to correct midfacial deformities presenting the clinical picture of unpleasant esthetic facial contour, facial asymmetries or malocclusion. The surgical treatment includes the separation of the maxilla into free segments to enable its repositioning in the desired, pre-surgically planned position. Regarding the fixation of the adjusted segments, the use of titanium mini-plates and screws is referred as the gold standard (Coskunses et al., 2015; He et al., 2015; Pan and Patil, 2014; Ueki et al., 2006).

Typically, the maxillary segment position is planned and primarily reconstructed with articulated dental models made from plaster casts before surgery. As osteotomies are conventionally based on two-dimensional (2D) lateral telerradiographies the precise intraoperative adjustment of the segments using surgical splints is often challenging. To overcome this issue alternative treatment approaches have been introduced. Preoperative virtual surgery planning and rapid prototyping surgical guides have been applied to ensure three-dimensional (3D) planning and separation of the segments in the exact position (He et al., 2015; Hirsch et al., 2009; Li et al., 2013; Mazzoni et al., 2015; Philippe, 2013). However, commercial straight titanium mini-plates, used for fixation still demand contouring to fit segmental maxillary geometry profiles for each individual patient, encountering a risk of inaccurate re-fixation of the segments (He et al., 2015). Furthermore, contouring of the titanium plates often comes along with repeated

bending leading to less stress resistance of the plate, increasing the risk of fatigue failure (Philippe, 2013). Custom made prefabricated titanium mini-plates have been investigated and discussed to allow precise control of the surgical procedure and decrease operative time (Mazzoni et al., 2015; Philippe, 2013). Beside titanium plates, poly-L-lactic acid plates and wires have been successfully used to achieve adequate postoperative maxillary stability (Egbert et al., 1995; Ueki et al., 2012). However, metal free solutions are not frequently used, titanium remains the material of choice though its removal is often indicated due to unclear potential bioactive corrosive products (Bianco et al., 1996; Stejskal and Stejskal, 1999; Weingart et al., 1994). Although zirconia or zirconium dioxide ceramic (ZrO₂) gained increasing popularity in the field of oral implantology its use as an osteosynthesis systems has not yet been described (Kubasiewicz-Ross et al., 2017; Manicone et al., 2007). The fact that ZrO₂ wear particles induced less pro-inflammatory gene expression compared to titanium particles is considered one rationale for its use (Obando-Pereda et al., 2014). Furthermore, ceramic-based miniplates are identified to cause less artifacts in magnetic resonance imaging (MRI) or computer tomography (CT) compared to titanium (Neumann et al., 2006). Although ZrO₂ is not resorbable it would fulfill the request for a bioinert and metal-free osteosynthesis system in maxillofacial surgery (Hayashi et al., 1992; Kubasiewicz-Ross et al., 2017).

The applicability of ZrO₂ osteosynthesis system to stabilize maxillary segments following a Le-Fort I advancement surgery has not yet

* Corresponding author.

E-mail address: l.hingsammer@gmail.com (L. Hingsammer).



Fig. 1. The experimental set up showing the two cast blocks, reflecting the separated maxillary segments stabilized in the desired position using virtually planned and individual fabricated ZrO₂ mini-plates.

Table 1
Mechanical properties of cortical bone, ZrO₂ and titanium in finite element analysis.

Material	Young's modulus (ε) GPa	Poisson ratio (ν)
Cortical bone	14.8	0.3
ZrO ₂	210	0.3
Titanium	144	0.33

Table 2
Number of elements and nodes for different structures and the overall model.

Part	No. of elements	No. of nodes
Bone	113,739	161,428
4 hole miniplate	5760	7884
3 hole miniplate	9725	13,152
Screw	4873	3132
Overall model	295,477	420,939




Fig. 2. Black arrows on the generated FEA model indicate the direction of loads (120 N, 500 N) simulating occlusal forces.

Table 3
Maximum Von Mises stress, P_{max} stress, yield strength and safety factor before plastic deformation under different loads for titanium and zirconium dioxide miniplates.

Load (N)	Material	P _{max} stress (MPa)	Max. Von Mises stress (MPa)	Yield strength (MPa)	Safety factor before plastic deformation
120	ZrO ₂	41.09	47.75	390	8.16
120	Titanium	35.18	47.87	828	17.29
500	ZrO ₂	72.83	182.82	390	2.13
500	Titanium	60.53	183.32	828	4.51

been evaluated. ZrO₂ belongs to the materials with the highest strengths suitable for medical use (von Wilmsowky et al., 2014). The aim of this study is to test if individual ZrO₂ mini-plates can stand occlusal forces and might constitute an appropriate solution for maxillary stabilization after Le Fort I osteotomy. Postoperative biomechanical behavior and stress distribution on titanium versus ZrO₂ mini-plates after Le Fort I advancement surgery was evaluated using 3D finite element analysis (FEA).

2. Materials and methods

Using two cast blocks, reflecting the separated segments of the maxilla, advancement of 5 mm and an extrusion of 2 mm was constructed. Individual fabricated ZrO₂ mini-plates and screws (Lithoz®, Vienna, Austria) were applied to stabilize the cast blocks in its position (Fig. 1). Two straight 3-hole and two angled 4-hole ZrO₂ mini-plates, secured with twelve cylindrical ZrO₂ screws all of which identical to the dimensions of the conventionally available titanium osteosynthesis system Modus 2.0 of Medartis®, Basel, Switzerland were used to achieve a stable fixation. Following this, 3D imaging of the cast blocks using a micro-CT computer tomography (RayScan 250E, Meersburg, Germany) with a voxel size of 65 μm was performed. The resulting DICOM data sets were then transferred to a stereolithography file format using VG Studio MAX 3.0 (Volume Graphics, Heidelberg, Germany). Abaqus CAE 6.12® software was used for creating 3D FE models of the ZrO₂ and titanium mini-plates, the cast blocks and the screws. All material properties including cortical bone were assumed to be isotropic, homogeneous and linear elastic. In Model I plates and screws were simulated to be made of conventional pure titanium grade IV and in model II of the existent printed ZrO₂. Material property assumptions regarding Young's modulus and Poisson ratios are listed in Table 1.




Fig. 3. (a) Three-dimensional highest Von Mises stress locations occurring in model I on the titanium screws (numbered from top to bottom) and plates. (b) Three-dimensional highest Von Mises stress locations on ZrO₂ screws and miniplates (model II) under 120 N of loading.

Values for ZrO₂ were obtained from Lithoz, Austria and that of pure titanium grade IV as well as that of bone were adopted from the literature (Ataç et al., 2008). Titanium has a pronounced yield behavior whereas ZrO₂ is a brittle material hardly allowing plastic deformation. Therefore, the maximum allowable value for Von Mises stress of ZrO₂ was set according to the yield strength value listed in Table 3. The number of elements and nodes identical for both models are listed in Table 2. The plate-to-screw, plate-to-bone and screw-to-bone interface assumed a full bonded condition to exclude micro-movements and to allow stress transfer continuity. The boundary conditions of screw – plate, plate-mono-cortical, and screw-mono-cortical fixation were created as hard contact surface condition. 120 N, 500 N were applied vertically in the molar and premolar region (Fig. 2).

Abaqus CAE 6.12® standard/implicit finite element solver served to analyze the FE models and Abaqus CAE 6.12® visualizer was used for visualizing the resulted stress concentrations.

3. Results

The 3D-FEA method was used to assess the Von Mises stress as well as the principal maximum stress (P_{max}) on zirconium dioxide (ZrO₂) and titanium miniplates in the Le Fort I advancement model. The maximum values of P_{max} and Von Mises stress of the miniplates under 120 N and 500 N vertical loading conditions are shown in Table 3 for model I (titanium) and II (ZrO₂).

All stress values are given in MPa (Newton per millimeter square). A color scale with 13 stress values served to evaluate quantitatively the stress distribution in the plates, screws and the adjacent bone tissue and to provide clear visualization of the stress concentrations.

3.1. Von Mises stress values in model I (titanium)

The resulted Von Mises stresses in the titanium miniplates under

500 N loading condition are shown in Fig. 3a. Maximum Von Mises Stress was 47.87 MPa (120 N) (Fig. 4a) and 183.32 MPa (500 N) respectively.

3.2. Von Mises stress values in model II (ZrO₂)

The simulation results indicate that the maximum von Mises stresses for the ZrO₂ miniplates were 47.75 MPa under 120 N and 182.82 MPa under 500 N vertical load.

Fig. 3b illustrates the location of Von Mises stresses of the ZrO₂ miniplates at 500 N.

In both models highest von Mises stresses were determined at the bending of the anterior placed plates. P_{max} for the ZrO₂ plates was 41.09 MPa under weak (120 N) and 72.83 MPa under heavy loads (500 N), respectively. Titanium miniplates resulted in P_{max} levels of 35.18 MPa for 120 N and 60.53 MPa for 500 N loaded models.

The safety factor before plastic deformation at titanium and crack at ZrO₂ occurs was 2.13 for ZrO₂ and 4.51 for titanium under a load of 500 N. When 120 N were applied, calculations revealed a safety factor of 8.16 for ZrO₂ and 17.29 for titanium (Figs. 4b and 5b).

4. Discussion

The present study evaluates the applicability of a 3D printed ZrO₂ osteosynthesis system to stabilize maxillary segments following a Le-Fort I advancement surgery. Results reveal a safety factor before plastic deformation respectively crack of 2.13 for ZrO₂ and 4.51 for titanium miniplates under a occlusal load of 500 N. Although no data regarding the bite force after Le Fort I Osteotomy exist, Harada et al. reported mean bite forces of 66.5 N after 2 weeks, 128.8 N after 4 weeks and 301.5 N after 6 months following BSSO surgery. Thus, occlusal loads of 120 N and 500 N were considered to simulate realistic bite forces during and after the bone healing phase of maxillary segments (Harada




Fig. 4. (a) Three-dimensional highest Von Mises stress locations occurring in model I on the titanium screws (numbered from top to bottom) and plates. (b) Three-dimensional highest Von Mises stress locations on ZrO₂ screws and plates (model II) under 500 N of loading.

et al., 2000). Multiple studies report various options, including titanium plates, poly-L-lactic acid plates and wires to achieve satisfactory maxillary stability after single-piece maxillary impactions and/or advancements. All methods provided satisfactory results without any appreciable differences (Egbert et al., 1995; Proffit et al., 1996, 1991; Skoczylas et al., 1988; Ueki et al., 2012). However, the 4-titaniumplate fixation technique constitute the gold standard as compared to the 2-plate fixation technique it significantly reduces stress on healing bones (Ataç et al., 2008). The usage of ZrO₂ devices has not yet been described.

Within the limitations of the study, a safety factor of 2.13 for the ZrO₂ miniplates under excessive occlusal load is considered adequate to ensure an uneventful clinical usage.

It is to be noted that to achieve comparable results, the evaluated ZrO₂ miniplates were fabricated according to the design and dimensions of conventionally available titanium miniplates used for model I without any pre-optimization. Thus, modifications of design and structure to even improve biomechanical behavior especially fracture resistance are clearly seen as new assignments for further investigations.



Titanium has been the gold standard for decades and still remains the material of choice for rigid fixation of freed maxillary segments (Coskunses et al., 2015; Philippe, 2013). Nevertheless, higher concentrations of titanium have been detected within tissue attached to implant surfaces and regional lymph nodes (Bianco et al., 1996; Weingart et al., 1994). These findings raise the question if titanium or its corrosive products have an impact on patients' individual health (Stejskal and Stejskal, 1999; Valentine-Thon and Schiwara, 2003). Although this issue is not fully elucidated, metal free solutions are not

routinely used. Zirconium ceramic has shown excellent biocompatibility and tissue integration (Manicone et al., 2007). Concerning dental implants the low affinity of zirconia to bacteria come along with adequate osseointegration of the material (Al-Radha et al., 2012; von Wilmsowsky et al., 2014).



Following the results of this study, ZrO₂ is considered as a potential solution to the raising claim of patients to stable, metal-free and bioinert osteosynthesis material without the routine need of removal (Pan and Patil, 2014; Verweij et al., 2016). Another crucial advantage of a ZrO₂ compared to a titanium osteosynthesis-system is that it does not cause artifacts in the CT or MRI scans (Neumann et al., 2006). This allows proper radiologic assessment of tissue adjacent to the osteosynthesis material and therefore unrestricted detection of diseases is possible (Neumann et al., 2006). Furthermore, patient specific fabricated osteosynthesis and the use of surgical cutting templates is considered to allow a precise positioning of the segments in the virtually planned position (He et al., 2015; Philippe, 2013).

Although clinical trials are missing, the findings of this study assume that individual ZrO₂ osteosynthesis screws and plates are suitable to stabilize freed maxillary segments after Le-Fort I advancement osteotomy. Further, a variety of possible indications for the use of ZrO₂ osteosynthesis material, especially in the field of facial traumatology, is considered. Nevertheless, ZrO₂ has gained increasing popularity as an implant material. Not only the fact that patients assert the claim to best esthetic results and with increasing frequency to metal-free solutions but also the bioinert property, especially when coated with saliva are reasons for its utilization.

However, increased fracture risk compared to titan implants due to low fracture toughness and stress shielding as a result of a very high



a



b

Fig. 5. Resulting Von Mises stresses illustrated during 500 N applied load. (a) titanium, (b) ZrO₂.

elastic modulus (210 GPa) compared to cancellous bone (14.8 GPa) are drawbacks of ceramics (von Wilmsowky et al., 2014). Furthermore, precise preoperative planning of the desired position of segments is essential as individual ZrO₂ implants would not allow any intraoperative adjustments compared to conventional titanium miniplates (Araújo et al., 2015; He et al., 2015; Hirsch et al., 2009).

As this study evaluates the biomechanical behavior of a ZrO₂ osteosynthesis system in the stabilization of the maxillary segments following Le-Fort I osteotomy for the first time, it has to be considered as a thought-provoking impulse to establish metal-free solutions. Following the results, this study can be considered as a pioneer work to introduce a new material (ZrO₂) for individually patient specific designed 3D printed implants in maxillofacial surgery. Beyond a doubt the not realistic simulation of the maxillary segments and the load only applied in one direction do not properly represent the dynamic loading during function. In addition, bone was modeled as linearly elastic and homogeneous even though bone in reality is anisotropic and inhomogeneous.

For clinical usage advance computational work is necessary to allow precise fabrication of individual ZrO₂ osteosynthesis parts and the surgical templates to ensure a predictable osteotomy of the segments. Nevertheless, the reported findings indicate a possible scope of patient

specific fabricated metal-free osteosynthesis systems.

Conflict of interest and source of funding

The authors do not have any financial interests, either directly or indirectly, in the products or information listed in the paper.

References

- Al-Radha, A.S.D., Dymock, D., Younes, C., O'Sullivan, D., 2012. Surface properties of titanium and zirconia dental implant materials and their effect on bacterial adhesion. *J. Dent.* 40, 146–153. <http://dx.doi.org/10.1016/j.jdent.2011.12.006>.
- Araújo, M.M., Lauria, A., Mendes, M.B.M., Claro, A.P.R.A., Claro, C.A., de, A., Moreira, R.W.F., 2015. Analysis of residual stress and hardness in regions of pre-manufactured and manual bends in fixation plates for maxillary advancement. *Oral Maxillofac. Surg.* 19, 369–373. <http://dx.doi.org/10.1007/s10006-015-0504-7>.
- Ataç, M.S., Erkmen, E., Yücel, E., Kurt, A., 2008. Comparison of biomechanical behaviour of maxilla following Le Fort I osteotomy with 2- versus 4-plate fixation using 3D-FEA. Part 1: advancement surgery. *Int. J. Oral Maxillofac. Surg.* 37, 1117–1124. <http://dx.doi.org/10.1016/j.ijom.2008.10.004>.
- Bianco, P.D., Ducheyne, P., Cuckler, J.M., 1996. Local accumulation of titanium released from a titanium implant in the absence of wear. *J. Biomed. Mater. Res.* 31, 227–234. [http://dx.doi.org/10.1002/\(SICI\)1097-4636\(199606\)31:2<227::AID-JBM9>3.0.CO;2-P10.1002/\(SICI\)1097-4636\(199606\)31:2<227::AID-JBM9>3.0.CO;2-P](http://dx.doi.org/10.1002/(SICI)1097-4636(199606)31:2<227::AID-JBM9>3.0.CO;2-P10.1002/(SICI)1097-4636(199606)31:2<227::AID-JBM9>3.0.CO;2-P).
- Coskunes, F.M., Kan, B., Mutlu, I., Cilasan, U., Celik, T., 2015. Evaluation of prebent

- miniplates in fixation of Le Fort I advancement osteotomy with the finite element method. *J. Cranio-Maxillofac. Surg.* 43, 1505–1510. <http://dx.doi.org/10.1016/j.jcms.2015.07.004>.
- Egbert, M., Hepworth, B., Myall, R., West, R., 1995. Stability of Le Fort I osteotomy with maxillary advancement: a comparison of combined wire fixation and rigid fixation. *J. Oral Maxillofac. Surg.* 53 (243–8–9).
- Harada, K., Watanabe, M., Ohkura, K., Enomoto, S., 2000. Measure of bite force and occlusal contact area before and after bilateral sagittal split ramus osteotomy of the mandible using a new pressure-sensitive device: a preliminary report. *J. Oral Maxillofac. Surg.* 58 (370–3–4).
- Hayashi, K., Matsuguchi, N., Uenoyama, K., Sugioka, Y., 1992. Re-evaluation of the biocompatibility of bioinert ceramics in vivo. *Biomaterials* 13, 195–200.
- He, W., Tian, K., Xie, X., Wang, X., Li, Y., Wang, X., Li, Z., 2015. Individualized surgical templates and titanium microplates for Le Fort I osteotomy by computer-aided design and computer-aided manufacturing. *J. Craniofac. Surg.* 26, 1877–1881. <http://dx.doi.org/10.1097/SCS.0000000000001938>.
- Hirsch, D.L., Garfein, E.S., Christensen, A.M., Weimer, K.A., Saddeh, P.B., Levine, J.P., 2009. Use of computer-aided design and computer-aided manufacturing to produce orthognathically ideal surgical outcomes: a paradigm shift in head and neck reconstruction. *J. Oral Maxillofac. Surg.* 67, 2115–2122. <http://dx.doi.org/10.1016/j.joms.2009.02.007>.
- Kubasiewicz-Ross, P., Dominiak, M., Gedrange, T., Botzenhart, U., 2017. Zirconium: the material of the future in modern implantology. *Adv. Clin. Exp. Med.* 26, 533–537. <http://dx.doi.org/10.17219/acem/63794>.
- Li, B., Zhang, L., Sun, H., Yuan, J., Shen, S.G.F., Wang, X., 2013. A novel method of computer aided orthognathic surgery using individual CAD/CAM templates: a combination of osteotomy and repositioning guides. *Br. J. Oral Maxillofac. Surg.* 51, e239–e244. <http://dx.doi.org/10.1016/j.bjoms.2013.03.007>.
- Manicone, P.F., Rossi Iommetti, P., Raffaelli, L., 2007. An overview of zirconia ceramics: basic properties and clinical applications. *J. Dent.* 35, 819–826. <http://dx.doi.org/10.1016/j.jdent.2007.07.008>.
- Mazzoni, S., Bianchi, A., Schiariti, G., Badiali, G., Marchetti, C., 2015. Computer-aided design and computer-aided manufacturing cutting guides and customized titanium plates are useful in upper maxilla waferless repositioning. *J. Oral Maxillofac. Surg.* 73, 701–707. <http://dx.doi.org/10.1016/j.joms.2014.10.028>.
- Neumann, A., Unkel, C., Werry, C., Herborn, C.U., Maier, H.R., Ragoss, C., Jahnke, K., 2006. Prototype of a silicon nitride ceramic-based miniplate osteofixation system for the midface. *Otolaryngol. Head. Neck Surg.* 134, 923–930. <http://dx.doi.org/10.1016/j.otohns.2006.01.022>.
- Obando-Pereda, G.A., Fischer, L., Stach-Machado, D.R., 2014. Titanium and zirconia particle-induced pro-inflammatory gene expression in cultured macrophages and osteolysis, inflammatory hyperalgesia and edema in vivo. *Life Sci.* 97, 96–106. <http://dx.doi.org/10.1016/j.lfs.2013.11.008>.
- Pan, Z., Patil, P.M., 2014. Titanium osteosynthesis hardware in maxillofacial trauma surgery: to remove or remain? A retrospective study. *Eur. J. Trauma Emerg. Surg.* 40, 587–591. <http://dx.doi.org/10.1007/s00068-013-0348-5>.
- Philippe, B., 2013. Custom-made prefabricated titanium miniplates in Le Fort I osteotomies: principles, procedure and clinical insights. *Int. J. Oral Maxillofac. Surg.* <http://dx.doi.org/10.1016/j.ijom.2012.12.013>.
- Proffit, W.R., Phillips, C., Prewitt, J.W., Turvey, T.A., 1991. Stability after surgical-orthodontic correction of skeletal Class III malocclusion. 2. Maxillary advancement. *Int. J. Adult Orthodon. Orthognath. Surg.* 6, 71–80.
- Proffit, W.R., Turvey, T.A., Phillips, C., 1996. Orthognathic surgery: a hierarchy of stability. *Int. J. Adult Orthodon. Orthognath. Surg.* 11, 191–204.
- Skoczylas, L.J., Ellis, E., Fonseca, R.J., Gallo, W.J., 1988. Stability of simultaneous maxillary intrusion and mandibular advancement: a comparison of rigid and nonrigid fixation techniques. *J. Oral Maxillofac. Surg.* 46, 1056–1064.
- Stejskal, J., Stejskal, V.D., 1999. The role of metals in autoimmunity and the link to neuroendocrinology. *Neuro Endocrinol. Lett.* 20, 351–364.
- Ueki, K., Marukawa, K., Shimada, M., Nakagawa, K., Alam, S., Yamamoto, E., 2006. Maxillary stability following Le Fort I osteotomy in combination with sagittal split ramus osteotomy and intraoral vertical ramus osteotomy: a comparative study between titanium miniplate and poly-L-lactic acid plate. *J. Oral Maxillofac. Surg.* 64, 74–80. <http://dx.doi.org/10.1016/j.joms.2005.09.015>.
- Ueki, K., Okabe, K., Moroi, A., Marukawa, K., Sotobori, M., Ishihara, Y., Nakagawa, K., 2012. Maxillary stability after Le Fort I osteotomy using three different plate systems. *Int. J. Oral Maxillofac. Surg.* 41, 942–948. <http://dx.doi.org/10.1016/j.ijom.2012.02.023>.
- Valentine-Thon, E., Schiwar, H.W., 2003. Validity of MELISA for metal sensitivity testing. *Neuro Endocrinol. Lett.* 24, 57–64.
- Verweij, J.P., Hassing, G.J.M., Fiocco, M., Houppermans, P.N.W.J., van Merkesteyn, J.P.R., 2016. Removal of osteosynthesis material because of symptoms after Le Fort I osteotomy: a retrospective study of 158 patients. *J. Cranio-Maxillofac. Surg.* 44, 1909–1912. <http://dx.doi.org/10.1016/j.jcms.2016.09.009>.
- von Wilmsow, C., Moest, T., Nkenke, E., Stelzle, F., Schlegel, K.A., 2014. Implants in bone: Part I. A current overview about tissue response, surface modifications and future perspectives. *Oral Maxillofac. Surg.* 18, 243–257. <http://dx.doi.org/10.1007/s10006-013-0398-1>.
- Weingart, D., Steinemann, S., Schilli, W., Strub, J.R., Hellerich, U., Assenmacher, J., Simpson, J., 1994. Titanium deposition in regional lymph nodes after insertion of titanium screw implants in maxillofacial region. *Int. J. Oral Maxillofac. Surg.* 23, 450–452.